



LARGE FRAGMENT IMPLANTS SYSTEM	
IMPLANT TRAY-1	
LC-LCP T. PLATE FOR Ø 4.5MM	
LC-LCP L. BUTTRESS PLATE FOR Ø 4.5MM- LEFT	
IMPLANT TRAY-2	
LC-LCP T. BUTTRESS PLATE FOR Ø 4.5MM	
LOCKING UPPER TIBIAL PLATE FOR Ø 4.5MM- LEFT	
IMPLANTS TRAY- 3	
LC-LCP DISTAL FEMORAL PLATE FOR Ø 4.5MM- LEFT	
IMPLANTS TRAY- 4	
LCP TROCHANTRIC PLATELEFT	
LCP MEDIAL DISTAL FEMUR PLATE (TOMOFIX)- LEFT	
IMPLANTS TRAY- 5	
LC-LCP T. PLATE FOR Ø 4.5MM	
LC-LCP L. BUTTRESS PLATE FOR Ø 4.5MM- RIGHT	
IMPLANTS TRAY- 6	
LC-LCP T. BUTTRESS PLATE FOR Ø 4.5MM	
LOCKING UPPER TIBIAL PLATE FOR Ø 4.5MM- RIGHT	
IMPLANTS TRAY- 7	
LC-LCP DISTAL FEMORAL PLATE FOR Ø 4.5MM- RIGHT	
IMPLANTS TRAY- 8	
LCP TROCHANTRIC PLATERIGHT	
LCP MEDIAL DISTAL FEMUR PLATE (TOMOFIX)- RIGHT	
SCREW TRAY	
CORTICAL SCREW (SELF TAPPING) Ø 4.5 MM	
LC-LCP SCREWS Ø 5.0 MM	
CANCELLOUS SCREWS Ø 6.5MM (F.TH)	
CANCELLOUS SCREWS Ø 6.5MM/16 MM TH	
CANCELLOUS SCREWS Ø 6.5MM/32 MM TH	
LC-LCP CANCELLOUS (F.TH) SCREW Ø 5.0MM	

LARGE FRAGMENT INSTRUMENT SET		
CAT NO.	NAME OF INSTRUMENTS	QTY
BOTTOM BOX INSTRUMENTS		
437.210	REDUCTION FORCEPS POINTED 210 MM	1
1617.22	TORQUE LIMITING SCREW DRIVER 4.5 MM (HEXAGONAL)	1
1617.023	HEXAGONAL SCREW DRIVER 4.5 MM	1
TRAY INSTRUMENTS		
1696.22	LCP COMBINED DRILL SLEEVE 3.2/4.3 MM	1
910.102	LCP DRILL SLEEVE/GUIDE 3.2 MM	2
910.128	LCP DRILL SLEEVE/GUIDE 4.3 MM	2
426.032	DRILL BIT 3.2 X 125 MM	2
427.143	DRILL BIT 4.3 X 225 MM	2
427.132	DRILL BIT 3.2 X 225 MM	2
241.201	K. WIRE 2.0 X 150 MM	5
421.045	DEPTH GAUGE 4.5 MM	1
	LARGE FRAGMENT INSTRUMENT BOX	1

Specification

- All implants are available in TT & SS.
- All instruments are available in SS.

Indication for use

- The Large Fragment Instrument Set is used for open reduction and internal fixation of large bone fractures, especially in the femur, tibia, and humerus. It is ideal for treating diaphyseal fractures, osteotomies, non-unions, and in cases requiring strong fixation with 4.5 mm or 5.0 mm screws and plates.

