


# NAIL EXTRACTOR INSTRUMENT SET



| NAIL EXTRACTOR INSTRUMENT SET |                                 |
|-------------------------------|---------------------------------|
| SR. NO.                       | NAME OF INSTRUMENTS             |
| 1                             | END CAP WRENCH 180MM SW11       |
| 2                             | WRENCH Ø5                       |
| 3                             | WRENCH Ø5                       |
| 4                             | WIRE CUTTER                     |
| 5                             | OPEN WRENCH SW11                |
| 6                             | HEX. SCREWDRIVER SW 4.5         |
| 7                             | HEX. SCREWDRIVER SW 5.0         |
| 8                             | LOCKING SCREW EXTRACTOR 250MM   |
| 9                             | HEX. SCREWDRIVER SW 2.5         |
| 10                            | HEX. SCREWDRIVER SW 3.5         |
| 11                            | HEX. SCREWDRIVER SW 4.0         |
| 12                            | QUICK COUPLING T-HANDLE         |
| 13                            | EXTRACTOR CHUCK UNIVERSAL LARGE |
| 14                            | EXTRACTOR CHUCK M6/M10          |
| 15                            | EXTRACTOR CHUCK SW4.0/M8        |
| 16                            | EXTRACTOR CHUCK SW4.5/M8        |
| 17                            | EXTRACTOR CHUCK M10*1           |
| 18                            | EXTRACTOR CHUCK M8*1            |
| 19                            | EXTRACTOR CHUCK M8              |
| 20                            | EXTRACTOR CHUCK UNIVERSAL SMALL |
| 21                            | EXTRACTOR CHUCK M11             |
| 22                            | EXTRACTOR CHUCK M10             |
| 23                            | SCREWDRIVER M8/M10              |
| 24                            | HAMMER M10*1                    |
| 25                            | UNIVERSAL JOINT M10/M10*1       |
| 26                            | ALUMINIUM BOX                   |

## Specification

- Instruments are available in SS.

## Indication for use

- The Nail Extractor Instrument Set is intended for use in the surgical removal of intramedullary nails from long bones such as the femur, tibia, and humerus. It is used during revision procedures, nonunion cases, infection management, or when implant removal is clinically required.