



PELVIC IMPLANTS SYSTEM
PLATES TRAY
PELVIC RECON PLATE LOW PROFILE Ø 3.5MM- STRAIGHT
PELVIC J PLATE Ø 3.5MM-LEFT
PELVIC J PLATE Ø 3.5MM-RIGHT
QUADRILATERAL SURFACE PLATE 3.5SHORT
QUADRILATERAL SURFACE PLATE 3.5STANDARD
QUADRILATERAL SURFACE PLATE 3.5-LONG
SPRING PLATE 3.5
INFRAPECTINEAL PLATE-SMALL - LEFT
INFRAPECTINEAL PLATE-LARGE - RIGHT
SUPRA PECTINEAL PLATE
3.5MM PUBIC SYMPHYSIS PLATE
PELVIC RECON PLATE 88° Ø 3.5MM CURVED
PELVIC RECON PLATE DOUBLE ANGLE
PELVIC RECON PLATE SINGLE ANGLE
PELVIC RECON PLATE 108° Ø 3.5MM CURVED
SCREWS TRAY
CANCELLOUS SCREWS Ø 4.0 MM (P.TH)
CANCELLOUS SCREWS Ø 4.0 MM (F.TH)
LARGE CANNULATED SCREWS Ø 7.0MM/16 MM TH
CORTICAL SCREW SELF TAPPING Ø 3.5 MM
CORTICAL SCREW SELF TAPPING Ø 4.5MM

PELVIC INSTRUMENTS SET		
CAT. NO.	NAME OF INSTRUMENTS	QTY
900.141	REDUCTION FORCEPS, OBLIQUE,1#	1
900.142	REDUCTION FORCEPS, OBLIQUE,2#	1
900.143	PELVIC REDUCTION FORCEPS WITH POINTS	1
900.144	BONE HOLDING FORCEP (SMALL)	1
900.145	SCREW DRIVER (HEX) SW2.5 LONG	1
900.146	BONE HOLDING FORCEP (LARGE)	1
900.147	PELVIC REDUCTION FORCEPS,1#	1
900.148	PELVIC REDUCTION FORCEPS,2#	1
900.149	ACETABULUM REDUCTION FORCEPS	1
900.150	LOCATING DRIVER	1
900.151	DRILL BIT 2.5 X 225MM	1
900.152	COMBINATION DRILL GUIDE	1
900.153	PLATE BENDING WRENCH	1
900.154	SCREW TAP 3.5 T-HANDLE	1
900.155	REDUCTION FORCEPS WITH 3BALL TIPS	1
900.156	SHARP HOOK	1
900.157	T-SHARP HOOK	1
900.158	PLATE BONDING WRENCH	1
900.159	UNIVERSAL SCREW DRIVER	1
900.160	HOHMANN RETRACTOR	1
900.161	DRILL GUIDE SLEEVE	1

Specification

- All implants are available in TT & SS.
- All instruments are available in SS.

Indication for use

- The Pelvic Instrument Set is used for open and minimally invasive fixation of complex pelvic and acetabular fractures. It supports procedures involving reduction, stabilization, and plating of the pelvic ring and acetabulum using specialized retractors, clamps, and reduction tools.