




**POSTERIOR CERVICAL IMPLANTS SYSTEM**

- LATERAL MASS / POSTERIOR CERVICAL POLY SCREWS 3.5 MM
- LATERAL MASS / POSTERIOR CERVICAL POLY SCREWS 4.0 MM
- POSTERIOR CERVICAL (LATERAL MASS) POLY SCREWS 3.5 MM (P.TH) FOR C1/C2 (TITANIUM)
- SPINAL ROD -Ø 3.5 MM
- TAPPERED SPINAL ROD - Ø 5.5 /3.5 MM (TITANIUM)
- CERVICAL LOCKING OCCIPITAL PLATE (TITANIUM)
- OCCIPATAL PLATE SCREW (TITANIUM)

| POSTERIOR CERVICAL INSTRUMENTS SET |                                                       |     |
|------------------------------------|-------------------------------------------------------|-----|
| CAT. NO.                           | NAME OF INSTRUMENTS                                   | QTY |
| 1632.006                           | SINGLE-AXIAL SCREW WRENCH                             | 1   |
| 1632.007                           | PRELOADED WRENCH                                      | 2   |
| 1632.004                           | POSTERIOR-CERVICAL HEX-WRENCH                         | 1   |
| 1632.008                           | DRILL                                                 | 2   |
| 1632.001                           | DRILL GUIDE                                           | 1   |
| 1632.024                           | GENERAL HANDLE                                        | 2   |
| 1632.016                           | CERVICAL POSITIONING-PIN SLEEVE                       | 1   |
| 1632.055                           | LATER MASS / POSTERIOR CERVICAL IMPLANTS BOX          | 1   |
| 1632.022                           | TEMPLATE (Ø3)                                         | 1   |
| 1632.015                           | POSTERIOR CERVICAL EXPANDING AWL                      | 1   |
| 1632.023                           | PROBE                                                 | 1   |
| 1632.025                           | POSITIONING PIN                                       | 1   |
| 1632.020                           | POSTERIOR CERVICAL DISTRACTOR                         | 1   |
| 1632.012                           | INTERNAL ROD BENDER                                   | 1   |
| 1632.021                           | POSTERIOR CERVICAL ROD BREAKER                        | 1   |
| 1632.019                           | POSTERIOR CERVICAL COMPRESSOR                         | 1   |
|                                    | POSTERIOR CERVICAL MECHANICAL ROD                     | 1   |
| 1632.003                           | COUNTER-TORQUE WRENCH                                 | 1   |
| 1632.009                           | POSTERIOR CERVICAL ROD HOLDER                         | 1   |
|                                    | LAMINA DETACHER                                       | 2   |
| 1632.026                           | HOOK DELIVER                                          | 1   |
| 1632.005                           | MULTI-AXIAL SCREW WRENCH (CERVICAL POLY SCREW DRIVER) | 2   |
|                                    | POSTERIOR CERVICAL HOOK HOLDER                        | 1   |
| 1632.009                           | POSTERIOR CERVICAL ROD HOLDER                         | 1   |
| 1632.014                           | POSTERIOR CERVICAL TAP                                | 1   |
| 1632.018                           | DEPTH GAUGE                                           | 1   |

## Specification

- All instruments are available in TT.
- All instruments are available in SS.

.

## Indication for use

- The Posterior Cervical System is used for spinal fixation procedures involving the posterior cervical spine. It is indicated for conditions like cervical fractures, instability, deformities, tumors, and degenerative disorders, enabling screw-rod fixation and fusion from the back of the neck.

