



V.A. HAND FIXATION IMPLANTS SYSTEM
PLATES TRAY
VARIABLE ANGLE LOCKING Y-PLATE, 3 HOLES HEAD
VARIABLE ANGLE LOCKING T- PLATE PLATES 3 H HEAD
VARIABLE ANGLE LOCKING CONDYLAR PLATES 2 H HEAD
VARIABLE ANGLE LOCKING STRAIGHT ADAPTION PLATE
VARIABLE ANGLE LOCKING STRAIGHT RECON PLATE
VARIABLE ANGLE LOCKING STRAIGHT ADAPTION FINGER PLATE
VARIABLE ANGLE LOCKING PHALANGEAL HEAD PLATE
VARIABLE ANGLE LOCKING BOX PLATE
VARIABLE ANGLE LOCKING STRUT PLATE
VARIABLE ANGLE LOCKING FIRST METACARPAL PLATE, LATERAL
VARIABLE ANGLE LOCKING FRACTURE PLATE
VARIABLE ANGLE LOCKING ROTATION CORRECTION PLATE
VARIABLE ANGLE LOCKING FIRST METACARPAL PLATE, DORSAL
VARIABLE ANGLE LOCKING METACARPAL NECK PLATE
VARIABLE ANGLE LOCKING ROTATION CONDYLAR PLATE
VARIABLE ANGLE LOCKING WEB PLATE
VARIABLE ANGLE LOCKING PHALANGEAL BASE PLATE
VARIABLE ANGLE LOCKING ROTATION T - PLATE
VARIABLE ANGLE LOCKING ROTATION Y-PLATE
SCREWS TRAY
V.A-LCP SCREW 1.3 MM
V.A-LCP SCREW 1.5 MM
V.A-LCP SCREW 2.0 MM
CORTICAL SCREW 1.3MM STAR HEAD
CORTICAL SCREW 1.5MM STAR HEAD
CORTICAL SCREW 2.0MM STAR HEAD

V.A. HAND FIXATION INSTRUMENTSET

CAT. NO.	NAME OF INSTRUMENTS	QTY
1623.003	BENDING IRON FOR 2.0MM FOR MAXILLOFACIAL PLATE / MINI BONE PLATE (PAIR)	1
1617.116	STAR SCREW DRIVER WITH SILICON HANDLE T6(1.5MM TIP)	2
1617.206	STAR TORQUE LIMITING SCREW DRIVER Ø T6 FOR 1.5/2.0MM	1
423.140	QUICK COUPLING ADAPTOR WITH T HANDLE	1
1597.125	SKIN HOOK SINGLE TWO PRONGE FOR HAND SYSTEM	2
1763.127	DOUBLE ENDED ELEVATOR FOR HANDFIXATION SYSTEM 3.5MM	1
1763.129	DOUBLE ENDED ELEVATOR FOR HAND FIXATION SYSTEM 5.0MM	1
1623.403	COUNTER SINK-CANNULATED S.S. T-HANDLE 2.0 MM	1
1598.014	DEPTH GOUGE FOR 1.5MM SCREWS	1
1723.215	NEEDLE BONE NIBBLER DOUBLE ACTION, 150 MM, ANGULAR FOR HANDFIXATION	1
1723.513	REDUCTION FORCEPS- POINTED- 130MM RATCHET LOCK FOR HAND FIXATION	1
429.213	DRILL TAP WITH COUPLING END Ø 1.3 MM	1
429.215	DRILL TAP WITH COUPLING END Ø 1.5 MM	1
429.220	DRILL TAP WITH COUPLING END Ø 2.0 MM	1
1617.306	STAR SCREW DRIVER WITH Q.C SHAFT- Ø T6 MM	1
1597.006	HOHMANN RETRACTOR, WIDTH 06MM	2
1596.251	LCP VARIABLE ANGLE DRILL SLEEVE 1.1MM (SINGLE SIDE) FOR 1.5MM SCREW	1
1596.255	LCP VARIABLE ANGLE DRILL SLEEVE 1.5MM (SINGLE SIDE)FOR 2.0MM SCREW	1
	DRILL BIT Ø 0.8MM	2
	DRILL BIT Ø 1.1MM	2
	DRILL BIT Ø 1.5MM	2
	ALUMINIUM BOX WITH TRAY	1

Specification

- All implants are available in TT & SS.
- All instruments are available in SS.

Indication for use

- The V.A. (Variable Angle) Hand Fixation Instrument Set is used for the fixation of fractures, fusions, and osteotomies in the hand and wrist. It allows precise placement of variable angle locking screws in small bone plates, offering stable fixation for complex hand trauma and reconstructive procedures.